



## pEco™ -BAD-nHis, PCR cloning Kit User Manual (Patent pending)

### Cloning PCR products for E Coli expression of N-term His-tagged protein

Cat#	Contents	Amounts	Application
IC-1003	pEco-BAD-nHis vector built-in Eco™ Cloning cells	10 tubes x 50ul/ea (for 10 rxn)	E Coli expression of N-term His-tag protein, best for toxic protein.
	Positive PCR insert	1 x 10ul/ea	
	20% L-arabinose	1ml	
	Sequencing primer pair	Forward and reverse 15ul/each, (25ng/ul)	

#### Storage:

Eco™ Cloning Kit is shipped on dry ice. Upon received, stored at -80°C. Once thawed, must be used, do not re-freeze. Product should be stable for 6 months.

#### Product Description:

##### 1. Introduction:

The revolutionary Eco Fusion *in vivo* cloning method is the easiest PCR cloning method available:


1. Simply amplify your gene of interest with a primer pair that is flanked with short arms homologous to the expression vector
2. Add 1µl of purified PCR into the engineered, vector build-in cloning cells
3. Immediately proceed to transformation.

##### 2. How it works:

The engineered *E Coli* strain in **GenTarget's Eco™ PCR Cloning Kit** has an enhanced *E Coli* competent cells enabling an *in vivo* joining reaction for cloning with no tube reactions.

#### Let the E Coli do the job for you *In Vivo!*

GenTarget provides *E Coli* cloning cells with a selection of built-in vectors for mammalian or *E Coli* expression systems. A proprietary process for making ready-to-use, *E Coli* cells with built-in vectors ensures low background and a positive cloning rate of greater than 90%.



**pEco-BAD-nHis** cloning cells has a built-in non-T7 based expression vector with an E Coli araBAD promoter. The promoter is tightly inhibited by an AraC regulatory protein expressed from the vector backbone. This basal level self-control mechanism is ideal for the expression of toxic proteins, avoiding the difficulties associated with direct cloning and expression. Target expression is induced by L-arabinose which releases the AraC regulatory protein, and the expression level can be optimized and protein solubility enhanced by varying the concentration of L-arabinose. The PCR insert is cloned in-frame with an N-terminal His-tag.

### 3. Key Features:

- The **easiest and most cost-effective** PCR cloning method available. Simply add 1 $\mu$ l of PCR insert into provided cells for transformation regardless of the insert's size and concentration
- No need to buy vectors and no tedious bench work preparing a vector backbone
- No need to buy cloning competent cells
- No need for any enzymes or any tube reactions
- Precise **directional cloning** of PCR products without any extra amino acids except the affinity tag (N-His)
- Flexibility to allow addition of any cleavage site for removal of N-terminal His-tag if desired
- Compatibility with any PCR product with or without a 3' -A overhang (the extra -A overhang, if it exists, will be removed in the cloning step)



- Can be used with PCR products of varying sizes, from 200 bp to 10 kb. The **same PCR product** can be used to construct multiple different expression vectors
- Engineered *E Coli* and mammalian expression vectors for high protein yields
- Great for high throughput cloning

#### 4. Protocol Outline:

Produce and clean PCR products



Add 1-2 µl of PCR product into the cloning cells provided, immediately proceed to transformation



Pick colonies, isolate plasmids by miniprep to verify the positive clones



Express protein from the saved glycerol stock

#### 5. Detailed Protocol:

##### 1. PCR primer design:

- ✿ PCR primers used for generating inserts for Eco<sup>TM</sup> Cloning must contain a 20 - 25bp homologous sequence corresponding to the built-in vector. Design your primer pair as follows:

**Fwd:** 5'- catggcgcacatcaccatcatcatcat + 20bp of (5'end gene specific forward sequence)

**Rev:** 5'- ttgtagcagggttaacacgcgtcta + 20bp of (3'-end gene specific reverse sequence)

- ✿ A protein cleavage site may be included in the forward primer to allow excision of the N-term tag if desired. Its codon sequences must be in frame and set between the homologous leader and the 20bp gene specific sequence.
- ✿ An **example** of PCR primer design:

To design the primer pair for the following gene sequence:



**atggcctctgtgaaggaaaatcc**actctagtcctacctgcatttctcagccttgct  
tacctgttgccaacattgggccaaccgaattcttccaatctttatcttggtgcca  
gcgagatgtcctcaacaaggagctgatgcagcagaatgggattggttatgtgta  
aatgccagcaata**cctgtccaagcctgactttta**

The PCR primer for vector **pEco-BAD-nHis** will be:

Fwd: 5'- catggcgcacaccatcatcatcat**atggcctctgtgaaggaaa**

Rev: 5'- ttgtagcaggtaacacgcgtcta**aaagtcaggcttggacagg**

If inserting a protein cleavage site, the Forward primer will be:

Fwd: 5' catgcatcatcaccatcatcatcatNNNNNN**gcctctgtgaaggaaaatcc**

(Where the NNNNNN is the in-frame codon sequence of the cleavage site).

### ☀ Notes:

1. GenTarget's cloning kits with the same terminal tags share PCR insert sites. The three Eco<sup>TM</sup> cloning kits with N-terminal tags, **Cat# IC-1001, IC-1002 and IC-1003**, can share the same PCR insert, and the two cloning kits with C-terminal tags, **Cat#: IC-1006 and IC-1007**, can share the same PCR insert.
2. A stop codon does not need to be included in the PCR reverse primer since a stop codon is already built in immediately after the PCR insert.

## 2. Target amplification by PCR:

- ☀ Amplify your target using any PCR amplification protocol that works for you. To minimize PCR errors, we recommend using high fidelity DNA polymerase.
- ☀ Use any PCR purification column to clean your PCR products. If you do not obtain a single, discrete band from PCR, gel purify your fragment.
- ☀ **Important:** if your PCR template can generate background clones having Amp resistance, treat the PCR product with DPNI or perform gel purification.

## 3. Transformation:



- ✿ Thaw Eco™ Cloning cells in ice-water.
- ✿ After they are completely thawed, add 1-2 µl purified PCR product (from 50ng to 150ng) into each vial of cells, and mix briefly by tapping the tube with your finger.
- ✿ For control vials, add 1µl positive PCR-insert (provided) as a positive control, and then add µl water to a negative control cells vial.
- ✿ Put tubes back on ice and proceed to heat shock at 42 °C for 40 seconds. (Note: Do not leave DNA-cells mixture on ice for prolonged period, less than 15min are fine). Put tubes back on ice for 1 min, add **250 µl** of SOC medium, and incubate at 37 °C, shaking for 1hr.
- ✿ **Plating:** take **250 µl** aliquot, spread out on pre-warmed LB-agar plates containing 100 µg/ml ampicillin. Grow colonies at 37 °C overnight.
- ✿ **Note:** In the absence of a PCR-insert, cells usually form background colonies; the no-insert negative control also generates a few colonies. In the presence of PCR insert, however, > 90% colonies are positive. Colony number varies depending on the quality and quantity of the PCR products. The concentration of purified PCR product can be from 20 ng/µl to 150 ng/µl with sizes ranging from 200 bp to 10 kb. For simplicity (and particularly for high throughput cloning) we recommend adding 1-2 µl of PCR product into the cloning cells. Regardless of the PCR product's concentration and size, it will generate enough colonies (5 ~ 100 colonies in general) for downstream work.

#### 4. **Save glycerol stocks for later expression and verification of positive clones:**

- ✿ Pick 2-5 colonies; propagate in LB/Amp, and incubate at 37 °C overnight
- ✿ Save an aliquot of each clone in LB-Glycerol medium containing 100 µg/ml ampicillin at a final concentration of 15% Glycerol.
- ✿ Isolate the plasmid DNAs using a DNA miniprep kit
- ✿ Confirm the positives by restriction digestion:
  - i. The PCR insert can be cut out at two unique sites: EcoRI + HpaI
  - ii. Run a 1.2% agarose gel. You should see two bands: 4 kb backbone + the PCR insert (or multiple bands when the cuts exist within the PCR-insert).
- ✿ Final sequencing verification. Use the provided sequencing primer pair. The sequencing primer comes in a ready-to-use dilution: use 1µl for each sequencing reaction with 500ng plasmid in 20µl volume.



Cat #	Vector	Forward primer	Reverse primer
IC-1003	pEco-BAD-nHis	IC-1003-fwd 5'- gaaataattttgttaact	IC-1003-rev 5'- gattaatctgtatcagg

## 5. Protein expression:

Once positive clones are confirmed, they can be used directly for protein expression without re-transformation into another strain.

- ✿ Add 10 µl of your positive clones in 4ml of LB medium with 100µg/ml ampicillin, grow at 37 °C overnight, shaking (225–250 rpm).
- ✿ The next day, measured OD should be approximately 1-2. Inoculate a large volume by making a 1:100 dilution of overnight culture in LB or SOB medium containing 100 µg/ml ampicillin. Grow the cultures at 37 °C with shaking to an OD600 = ~0.5.
- ✿ Induce expression by adding L-arabinose in a range of final concentrations (from 0.2% to 0.00005%) to find the concentration for optimal expression.
- ✿ Remove a 1 ml aliquot of cells from each tube for analysis of protein expression. Be sure to save aliquots of uninduced control samples.
- ✿ Analyze protein expression by SDS-PAGE or other methods.
- ✿ **Harvest** cells by centrifugation.
- ✿ Lyse the cell pellet using lysis reagent. Following lysis, run a gel protein analysis.
- ✿ **Purification:** use your favorite protocol and reagent to purify the expressed His-tagged protein by affinity chromatography.

## 6. Vector maps:

The figure below summarizes the vector map of pEco-BAD-nHis. The **complete nucleotide sequence is** available for download from our Support page (<http://www.gentarget.com/support/vector-sequences/>). To make your clone map, simply paste your gene sequence (**not included the flanking sequences of both ends**) in the Red highlighted position (replacing the **NNNN..NN**). In most case, the pasted sequence is: "ATG...to...last codon".




## Cloning site for pEco-BAD-nHis vector

```

                                EcoRI
                                ~~~~~
301  CTAGAAATAA TTTTGTTTAA CTTTAAGAAG GAGAATTCAC CATGGCG CAT
      GATCTTTATT AAAACAAATT GAAATTCTTC CTCTTAAGTG GTACCGC GTA
                                HpaI
                                ~~~~~
      His-tag          PCR insert
351  CAC CAT CAT CAT CAT NNNNN NNNNTAGACG CGTGTTAACC TGCTAACAAG
      GTG GTA GTA GTA GTNNNNN NNNNATCTGC GCACAATTGG ACGATTGTTC
  
```

## Map of pEco-BAD-nHis vector



O2 operator: 4-19  
 CAP binding with 1/1/2: 161-251  
 Arabinose promoter: 248-276  
 RBS: 329-334  
 GOI cloning ends: 365-366  
 Amp: 901-1761  
 pUC19 pro: 1906-2579  
 araC gene: 3988-3110

## 7. Troubleshooting:

Problems	Solution
No colony	<ul style="list-style-type: none"> <li>☠ Be sure to set up a positive control transformation using the provided positive PCR insert1, which should give you 10~100 colonies.</li> <li>☠ Spread all of the transformation mixture onto the plate.</li> </ul>



Background colonies	<ul style="list-style-type: none"> <li>☀ Be sure to set up a background control plate in which no PCR product was added to the cells. It should generate 0 ~ 5 colonies or less than 10% compared to plates with the insert (Note: in the absence of a PCR insert, cells force vector self-ligation resulting in a few background colonies).</li> <li>☀ Make sure that the PCR's template does not cause background colonies; If it does, clean PCR products by gel-isolation or treatment with DPNI.</li> <li>☀ Plate less transformation mixture onto the plate.</li> </ul>
Satellite colonies	<ul style="list-style-type: none"> <li>☀ Be sure to use the right amount of antibiotics in the LB plate, and make fresh LB plates if necessary.</li> <li>☀ Use carbenicillin instead of ampicillin if applicable.</li> <li>☀ Do not incubate plates longer than 16 hours.</li> <li>☀ Try to avoid picking the tiny satellite colonies.</li> </ul>

## Related Products:

Cat#	Product Name	Amount	Application
<a href="#">IC-1001</a>	PCR cloning kit	kit	PCR cloning kit with a built-in vector (T7 promoter based) in provided cloning cells for <b>E Coli expression of N-term His-tagged protein.</b>
<a href="#">IC-1002</a>	PCR cloning kit	kit	PCR cloning kit with a built-in mammalian expression vector (with neomycin selection marker) in provided cloning cells. The vector containing an engineered super CMV promoter for high-yield <b>mammalian expression of N-term His tagged protein</b>
<a href="#">IC-1004</a>	PCR cloning kit	kit	PCR cloning kit with a built-in vector (T7 promoter based) in provided cloning cells for <b>E Coli expression of N-term GST-tagged protein.</b>
<a href="#">IC-1006</a>	PCR cloning kit	kit	PCR cloning kit with a built-in vector (T7 promoter based) in provided cloning cells, for <b>E Coli expression of C-term His-tagged protein.</b>
<a href="#">IC-1007</a>	PCR cloning kit	kit	PCR cloning kit with a built-in mammalian expression vector (with Neomycin selection marker) in provided cloning cells, for <b>mammalian expression of C-term His-tagged protein.</b>

## References:

- Oliner et al., 1993, Nucleic Acids Res. 1:5192-97
- Aslanidis et al., 1994, Genome Res. 4 :172-177
- Kaluz et al. Nucl. Acids Res..1992; 20: 4369-4370