



## FAQ for pre-made Adenovirus:

### 1. How do I use the GenTarget's pre-made Adenovirus?

Pre-made Adenovirus is provided ready to use. Simply add it into the mammalian cell line of your choice. The amount of virus to add is dependent upon your cell types, or your titer the virus amount based on MOI number.

### 2. What is MOI?

To obtain the desired expression level, you need a suitable multiplicity of infection (MOI). The MOI is defined as the number of virus particles per cell, and a range of MOIs, from 0, to 1000 should be tested to determine the MOI to use. Normally, actively dividing cells require less virus (a lower MOI) than non-dividing cells. An MOI that is too high may lead to cell death due to the toxicity of the virus or of the expressed target; an MOI that is too low may result in a low positive transduction rate and low target expression levels. You should use the minimal MOI required to produce the desired expression level or positive transduction rate (usually 100%). For most cell lines an MOI of 1-10 is fine.

### 3. What is the control virus?

The **AVP-Null** control virus (sold as a separate product) serves as the negative control for adenovirus treatment. The Null-control adenovirus is produced from an Adenovector cloned with a Null sequence which does not express any target.

### 4. How do I know the Adenovirus is working in my cells?

Adenovirus can transduce the majority of mammalian cells very well -- including human, mouse, and other species--whether dividing or non-dividing. However, there are a few cell types that cannot be transduced or that can be transduced only with very low efficiency. For those types of cells, you may wish to try lentivirus. The best way to find out if your cells are susceptible to adenoviral transduction is to use a fluorescent control adenovirus, which will let you easily visualize transduction as fluorescence signal. GenTarget provides pre-made adenoviruses, expressing a wide selection of fluorescent markers including GFP, RFP, YFP, and CFP, as well as fluorescent fusion markers.

### 5. How long will target expression last?

Typically, target expression can be detected after 12-24 hours. Depending upon the cell types or dividing cycles, the expression peaks at 2-3 days post



GenTarget Inc

7930 Arjons Drive, Suite B  
San Diego, CA 92126  
Phone: (858) 6788683  
Email: [orders@gentarget.com](mailto:orders@gentarget.com)  
Fax: (800) 3804198

-transduction and persists for as long as the viral genome is present, from 1 week to 6 months or longer.

#### 6. Can pre-made adenovirus be used *in vivo*?

Yes GenTarget provides adenovirus in PBS with 5% sucrose for *in vivo* applications. **NOTE:** all GenTarget pre-made Adenoviruses are **for research use only and not for clinical or therapeutic applications.**

#### 7. What buffer is the virus provided in?

GenTarget's adenoviruses are provided as 200 µl aliquots in two formats:

- Crude viral stock in DMEM medium with 10% serum
- Concentrated adenovirus re-suspended into PBS with 5% sucrose for *in vivo* use.

#### 8. What about bio-safety? What is RCA?

Adenoviruses are safe to use. They are non-replicable viral particles produced from an expression vector derived from human adenovirus type 5 genome in which the entire E1 (4.3kb) and E3 region was removed. E1 protein is essential for the viral replication without it, the packaged adenovirus cannot be replicated in target cells. The presence of the viral genome is transient since it will eventually be diluted out as cell division occurs. (Note: In E1 complementing cell lines, like 293A cells, the adenovirus can be replicated or amplified.)

Despite the adopted safety features, recombinant human adenovirus has been classified as a biosafety level II agent, and you will need a BL-2 level facility to work with it. It should be noted that cell culture facilities in most institutions are certified as BL-2 level. Please refer to CDC and NIH guidelines for more details about adenovirus handling.

In rare instances, the E1 sequence from the genome of the virus producing (293A) cells is integrated into and adenoviral trans-gene construct by sequence homologous recombination, resulting in a (**RCA**). In theory, this would most likely occur in large-scale virus amplification. Wild type, replication competent adenoviruses could cause cold symptoms but generally do not cause serious illness. All GenTarget adenoviruses are made on a very small scale and RCA testing is not routinely performed for pre-made adenovirus. When performing large scale adenovirus applications, we use PCR to measure the level of replication competent adenovirus. **Note:** When RCA occurs, it will quickly overtake the non-replicable virus and cause



**GenTarget Inc**

7930 Arjons Drive, Suite B  
San Diego, CA 92126  
Phone: (858) 6788683  
Email: [orders@gentarget.com](mailto:orders@gentarget.com)  
Fax: (800) 3804198

cell death. To avoid the occurrence of RCA, viruses should be produced and amplified in low passage 293A cells.

### **9. How is viral titer measured??**

The titer (IFU/ml) of GenTarget's adenovirus is measured via a fluorescent marker (GFP or RFP) after transduction of the virus into HeLa cells. This titer approximates the real infectious units measured by biological (plaque) assay, and is used as the reference titer for the non-fluorescent construct.

### **10. How are the adenovirus vectors constructed?**

GenTarget uses its proprietary Eco cloning technology (vector built-in cloning competent cells) to directly insert a target into the E1E3 deleted human adenovirus 5 genome.

### **11. Are there any antibiotic markers included in the adenovirus constructs?**

Adenovirus is not integrated into the host cell genome therefore, it is not stable for long term expression and we do not include an antibiotic marker in adenovirus constructs. In some cases, we do include a fluorescent marker under the same promoter for the specific target, mediated by F2A element. The fluorescent marker enables monitoring of viral performance and target expression, as well as selection via fluorescent signal. For long term stable expression, please use our pre-made lentivirus which features a wide variety of selection markers.

### **12. Adenovirus vs. Lentivirus**

Both Adenovirus and lentivirus can be transduced into dividing and non-dividing cells, and both are used as expression delivery tools for mammalian cell lines or primary cells. Unlike lentivirus, adenovirus is non-insertional so constructs do not integrate into the host genome and will not affect the activity of host genes. Adenovirus is a transient expression delivery vehicle; once transduced into mammalian cells, the recombinant adenoviral transgene target will be expressed until the viral genome is diluted by cell division. Transgene expression decreases gradually in actively dividing cells (normally in 1-2 weeks) but expression can persist for a longer time in non-dividing or slowly dividing cells, such as skeletal muscle cells or neurons. By contrast, lentivirus delivers stable, long term expression. Lentivirus also is less immunogenic to human cells and less toxic. Adenovirus may, however, have better transduction efficiency in some cell types.



**GenTarget Inc**

7930 Arjons Drive, Suite B  
San Diego, CA 92126  
Phone: (858) 6788683  
Email: [orders@gentarget.com](mailto:orders@gentarget.com)  
Fax: (800) 3804198

**Attachment: GenTarget's Pre-made Adenovirus Product categories:**

<b>Product Categories</b>	<b>Product Description (Click to see the products under the category)</b>
<b>Common Markers &amp; Enzymes</b>	<a href="#">Express a enzyme or maker, such as CRISPR /Cas9, Luciferase /GFP, or Luciferase /RFP,...</a>
<b>Fluorescent Controls</b>	<a href="#">Express enhance GFP, RFP, CFP, BFP and Null-control adenovirus....</a>
<b>iPS Stem Cell Factors</b>	<a href="#">Express the iPSC factors (OCT3/4, SOX2, LIN28)...</a>