



Lin28 & RFP stable cell line manual

Catalog Number	Amount	Storage
SC016	1 vial of cells (2 x 10 ⁶ cells) in 80% DMEM, 10% FBS, 10% DMSO	Liquid nitrogen

Product Description

The 293 Cell Line is a permanent line established from primary embryonic human kidney transformed with sheared human adenovirus type 5 DNA. The expressed E1A adenovirus gene allows these cells to produce very high levels of protein. Lin28 & RFP stable cell is transformed from the 293 cell line and stably expresses human Lin28 gene. It is established by transduction of Ln28 expression lentivirus. Lin28 and RFP are bicistronically expressed (mediated by Ferin 2A self-cleavage element) in high-level under [suCMV promoter](#). Lin28 is expressed with C-terminal His tag, and RFP has no any tags. Blasticidin marker is expressed under PGK promoter for selecting the stable cells.

Culture procedures

1. Thaw the frozen vial of cells quickly in a 37°C water bath (1~3min), decontaminate the outside of the vial with 70% ethanol.
2. Transfer the entire contents of the cryovial into a T-75 cm² flask containing 15 ml of pre-warmed complete medium. Incubate the cells overnight in a 37°C incubator, 5% CO₂.
3. The following day, replace the medium with 15 ml of pre-warmed, complete medium (Optional: add 10 µg/ml blasticidin in medium).
4. Incubate the cells and monitor cell density.
5. Pass cells (1:10 dilution) when the culture reaches 80-90% confluent.
6. Freeze cells at a density of 3 x 10⁶ cells/ml using 90% complete medium with 10% DMSO.

Complete medium

D-MEM (high glucose)
2mM L-glutamine
10% Fetal Bovine Serum (FBS)
0.1 mM MEM Non-Essential Amino Acids (NEAA)
1% Pen-strep (optinal)

Quality Control

Each vial contains greater than 2 x 10⁶ cells with >95% viability before freeze. Cells are tested free of bacteria, viruses, mycoplasma.

Warranty

This product is warranted to perform as described when used in accordance with this manual. Gentarget's sole remedy for breach of warranty should be, at the option of Gentarget, to repair or replace the product. **This product is for research use only.**